



3zahn Conical Connection
— Catalog —

Length
Diameter

8mm

10mm

11.5mm

13mm



| | | | |
|------------|-------------|--------------|-------------|
| ZDFCIP8375 | ZDFCIP10375 | ZDFCIP115375 | ZDFCIP13375 |
|------------|-------------|--------------|-------------|



| | | | |
|-----------|------------|-------------|------------|
| ZDFCIP842 | ZDFCIP1042 | ZDFCIP11542 | ZDFCIP1342 |
|-----------|------------|-------------|------------|




| | | | |
|-----------|------------|-------------|------------|
| ZDFCIP847 | ZDFCIP1047 | ZDFCIP11547 | ZDFCIP1347 |
|-----------|------------|-------------|------------|



| | | | |
|-----------|------------|-------------|------------|
| ZDFCIP852 | ZDFCIP1052 | ZDFCIP11552 | ZDFCIP1352 |
|-----------|------------|-------------|------------|





| | Height mm | Ref. |
|---|-----------|---------------|
|  | 3mm | ZAGM-203-3-RC |
|  | 4mm | ZAGM-203-4-RC |
|  | 5mm | ZAGM-203-5-RC |
|  | 6mm | ZAGM-203-6-RC |
|  | 7mm | ZAGM-203-7-RC |

Healing Caps - Ø4.5

Additional Info:

Recommended torque - manual only light finger - 5-10 Ncm. using hand driver 1.2mm.

Material: Titanium

| | Description | Ref. |
|---|-------------|-------------|
|  | Open Tray | ZAGM-301-RC |
|  | Close Tray | ZAGM-305-RC |

Implant Transfer - Ø4.5

Additional Info:

Recommended torque - manual only light finger tight (5-10 Ncm) using hand driver 1.2mm.

Material: Stainless Steel

| | Height mm | Ref. |
|---|-----------|-------------|
|  | 13mm | ZAGM-214-RC |

Implant analog

Material: Stainless Steel



| | Height mm | Ref. |
|--|-----------|----------------|
|  | 3mm | ZAGM-203-3W-RC |
|  | 4mm | ZAGM-203-W4-RC |
|  | 5mm | ZAGM-203-W5-RC |
|  | 6mm | ZAGM-203-W6-RC |
|  | 7mm | ZAGM-203-W7-RC |

Wide Healing Caps - Ø6.2

Additional Info:

Recommended torque - manual only light finger - 5-10 Ncm. using hand driver 1.2mm.

Material: Titanium

| | Description | Ref. |
|--|-------------|-------------|
|  | Hex | ZAGM-903-RC |
|  | Non-Hex | ZAGM-904-RC |

Temporary Abutment - Ø4.0

Additional Info:

Recommended torque - 20 Ncm.

Material: PEEK





| | Height mm | Ref. |
|--|-----------|-------------|
|  | 13mm | ZAGM-202-RC |




Direct Prosthetic Screw

Additional Info:

Recommended torque - 25 Ncm.

Material: Titanium

| | Height mm | Ref. |
|---|-----------|--------------|
|  | 1mm | ZAGM-101-1RC |
|  | 2mm | ZAGM-101-2RC |
|  | 3mm | ZAGM-101-3RC |
|  | 4mm | ZAGM-101-4RC |

| | Height mm | Ref. |
|---|-----------|--------------|
|  | 1mm | ZAGM-102-1RC |
|  | 2mm | ZAGM-102-2RC |
|  | 3mm | ZAGM-102-3RC |

Concave Shoulder Abutment - Ø4.5

Additional Info:

Recommended torque - 25 Ncm.




Material: Titanium

15° Esthetic Angulated Abutment - Ø4.5

Additional Info:

Recommended torque - 25 Ncm.

Material: Titanium


| | Height mm | Ref. |
|---|-----------|----------------|
|  | 1mm | ZAGM-103-E1-RC |
|  | 2mm | ZAGM-103-E2-RC |
|  | 3mm | ZAGM-103-E3-RC |






25° Esthetic Angulated Abutment - Ø4.5

Additional Info:

Recommended torque - 25 Ncm.

Material: Titanium



| | Height mm | Ref. |
|---|-----------|---------------|
|  | 1mm | ZAGM-700-1-RC |
|  | 2mm | ZAGM-700-2-RC |
|  | 3mm | ZAGM-700-3-RC |
|  | 4mm | ZAGM-700-4-RC |
|  | 5mm | ZAGM-700-5-RC |



Straight Multi Unit for Internal Hex - Ø5.0

Additional Info:

Recommended torque - 25 Ncm.

Material: Titanium

| | Height mm | Ref. |
|---|-----------|---------------|
|  | 1mm | ZAGM-717-1-RC |
|  | 3mm | ZAGM-717-3-RC |





17° Multi Unit Angulated - Ø5.0

Additional Info:

Recommended torque - 25 Ncm.

Material: Titanium

| | Height mm | Ref. |
|---|-----------|---------------|
|  | 1mm | ZAGM-730-1-RC |
|  | 3mm | ZAGM-730-3-RC |




30° Multi Unit Angulated - Ø5.0

Additional Info:

Recommended torque - 25 Ncm.

Material: Titanium

| | Description | Height mm | Ref. |
|---|--------------------|-----------|--------|
|  | Scan post for M.U. | 7mm | ZSBMU |
|  | Ti-Base | 0.3mm | ZTBMU |
|  | Digital Analog | 10mm | ZAD-MU |




Multi-Unit CAD/CAM

Additional Info:

Recommended torque - 25 Ncm.

Material: Peek , Stainless Steel






| | Description | Height mm | Ref. |
|---|-------------------------|-----------|-------------|
|  | M.U. transfer open tray | 22mm | ZAGM-508 |
|  | M.U. analog | 12mm | ZAGM-507 |
|  | M.U. healing cap | 5mm | ZAGM-509 |
|  | M.U. temporary cylinder | 12mm | ZAGM-505 |
|  | M.U. plastic cylinder | 12mm | ZAGM-506 |
|  | M.U. Screw | 3.8mm | ZAGM-506-SC |
|  | M.U. temporary cylinder | 12mm | ZAGM-512 |

Multi-Unit Prosthetics

Additional Info:

Recommended torque - 25 Ncm.

Material: Titanium, Stainless Steel
Plastic,PEEK

| | Height mm | Ref. |
|---|-----------|------------|
|  | 0.5mm | ZTBH-05-RC |
|  | 1.5mm | ZTBH-15-RC |
|  | 3mm | ZTBH-30-RC |

| | Height mm | Ref. |
|---|-----------|-------------|
|  | 0.5mm | ZTBNH-05-RC |
|  | 1.5mm | ZTBNH-15-RC |
|  | 3mm | ZTBNH-30-RC |

Engaging Ti-Base - Ø4.5

Additional Info:

Recommended torque - 25 Ncm.

Material: Titanium

Non-Engaging Ti-Base - Ø4.5

Additional Info:


Recommended torque - 25 Ncm.

Material: Titanium

|  | Height mm | Ref. |
|---|-----------|----------|
| | 10mm | ZSBRC-10 |

Scan post

Material: PEEK


|  | Height mm | Ref. |
|--|-----------|--------|
| | 9mm | ZAD-RC |

Digital Analog

Material: Stainless Steel



| | Height mm | Ref. |
|---|-----------|---------------|
|  | 1mm | ZAGM-106-1-RC |
|  | 2mm | ZAGM-106-2-RC |
|  | 3mm | ZAGM-106-3-RC |
|  | 4mm | ZAGM-106-4-RC |
|  | 5mm | ZAGM-106-5-RC |



Lock Attachment

Additional Info:
 Recommended torque - 25 Ncm.-
 Using hex driver 1.2mm.

Material: Titanium

Lock Set For RC Connection

| | | |
|---|--------------------------|-------------|
|  | Lock Plastic Cap-RC | ZAGM-106-P |
|  | Lock Soft Plastic Cap-RC | ZAGM-106-PS |
|  | Lock Protective Disk-RC | ZAGM-106-D |
|  | Lock Hard Plastic Cap-RC | ZAGM-106-PH |
|  | Lock Metal Housing-RC | ZAGM-106-M |

*Lock can be sold separately.






3zahn Monoblock implant offers the perfect solution for a unique one-piece implant with its flexible neck design that reaches **up to 25° angulation**.



The implant's versatility offers immediate loading in **multiple-unit restorations**.



Conveniently it can be placed in the tightest and **narrowest ridges**.




A wide range of diameters from 3.3mm up to the 4.7mm

Body Characteristics



Flexible Neck



Sharp and Deep Threads

Material: Titanium alloy Ti 6Al 4V ELI

Connection: One-Piece implant with an integrated abutment carrier

Item includes: Implant carrier

Surface: Sand blasting and Acid-Etching

| Diameter | Length | 6mm | 8mm | 10mm | 11.5mm | 13mm | 16mm |
|----------|--------|----------|----------|-----------|------------|-----------|-----------|
| Ø3.3 | | - | ZPTR833 | ZPTR1033 | ZPTR11533 | ZPTR1333 | ZPTR1633 |
| Ø3.75 | | ZPTR6375 | ZPTR8375 | ZPTR10375 | ZPTR115375 | ZPTR13375 | ZPTR16375 |
| Ø4.2 | | ZPTR642 | ZPTR842 | ZPTR1042 | ZPTR11542 | ZPTR1342 | ZPTR1642 |
| Ø4.7 | | ZPTR647 | ZPTR847 | ZPTR1047 | ZPTR11547 | ZPTR1347 | ZPTR1647 |

One-Piece:



One-Piece (Long)
ZITSL



One-Piece (Short)
ZITSLS

3Zahn Kit

Kit includes:



Torque wrench
ZTW

Implant Insertion Tools:




Handpiece 2.5mm (Long)
ZHAL



Handpiece 2.5mm (Short)
ZHAS



2.5mm (Long)
ZIT2.5



2.5mm (Short)
ZITS2.5


Prosthetic Drivers Hex 1.2mm:



Long
ZAGM-400



Short
ZAGM-400-S



Hand
ZAGM-400-R



Handpiece
ZHB

Tools:



Drill extender
ZDE



Parallel pin x3
ZPP
Height mm

Drill Stoppers:

Height mm

6mm 8mm 10mm 11.5mm 13mm

Drills:



ZDM



ZDB2.0



ZDB2.8



ZDB3.2




ZDB3.7




ZDB4.2



ZDB4.7



*Instruments can be order separately



3Zahn Kit

Ref. ZSKT



FDA CE

www.3zahn.com